
PROVOX[®] Dilator

CE 0413

Rx ONLY









ATOS



MEDICAL

Leading the way in ENT






Provox Dilator 17/20

REF 7122, 7123

2 / 250 ml (8,5 fl oz)

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden.

Contents

ENGLISH	6
DEUTSCH	15
FRANÇAIS	26
ITALIANO	36
ESPAÑOL	46
PORTUGUÊS.....	56
EESTI	66
MAGYAR	75
POLSKI.....	85
РУССКИЙ	95
日本語版.....	106
한국어	116
繁體中文	125
ORDERING INFORMATION.....	132

The manual, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

Intended use

The Provox voice prosthesis system is intended for use of prosthetic voice restoration after total laryngectomy.

The Provox Dilator 17 and the Provox Dilator 20 are used for dilating tracheo-esophageal (TE) punctures that shrink very fast or are too narrow to insert the selected Provox voice prosthesis. The dilator may also be used to temporarily block or stent the TE puncture.

For Home Care Use the dilator is intended for single patient re-use, and may only be used as part of the Provox NID system. It shall only be used by patients with sufficient manual dexterity, acceptable vision and satisfactory cognitive ability.

For Professional Use the dilator may be re-used by a clinician in clinic and shall be re-sterilized between patients.

Contraindications

- The device is not intended to be used by patients that have not received adequate training by their clinician.
- The device is not intended to be used for puncture dilation at the time of surgical creation of the puncture.
- Do not use the device in case of a small tracheostoma where the dilator could obstruct breathing.

Device description

The Provox Dilator 17 and the Provox Dilator 20 are tapered curved silicone rods used for dilating (increase the diameter of) TE punctures. The Provox Dilator 17 shall be used with Provox voice prostheses with an outer diameter of 17 Fr and the Provox Dilator 20 with outer diameter 20 Fr. The dilator shall only be used and be prescribed for patient use, by clinicians trained in the care and rehabilitation of laryngectomized patients.

 This product and its packaging are latex-free.

The Provox Dilator Figure 1

- 1.1: Medallion
- 1.2: Safety strap
- 1.3: Device handle
- 1.4: Notch
- 1.5: Retention collar
- 1.6: Tip

The size is marked on the medallion.

Dimensions:

Provox Dilator 17: The diameter is 11.5 Fr at the tip and increases to 19 Fr.

Provox Dilator 20: The diameter is 13 Fr at the tip and increases to 22 Fr.

Length between tip and notch: 60 mm/2.5 inches; length between retention collar and notch: 20 mm

Recommended Lubricants

The following water-soluble lubricants are recommended for use with the dilator:

KY Jelly®

SurgiLube®

Recommended anaesthetizing gel

The following anaesthetizing gel is recommended for use with the dilator:

Xylocaine®/(Lidocaine) Gel.

CAUTION: Lubricants or anaesthetizing gels that are not water soluble shall not be used at any time.

WARNINGS

General:

- Do not use the device if it shows any sign of structural damage. A damaged device may break and parts of it can enter the airways and have to be removed immediately.
- Do not use lubricating (anaesthetizing) gel if the patient has allergies related to these substances.
- Never remove the medallion from the device. The medallion is of utmost importance for the safety. It helps to reduce the risk of the device accidentally falling into the trachea (windpipe).

Professional use:

- To avoid cross contamination, use appropriate hygienic procedures and make sure that each used device is cleaned, disinfected, dried and steam sterilized.
- The material in the Provox Dilator withstands 50 cleaning and steam sterilization procedures at 137°C/278°F for 3 minutes or at 134°C/273°F for 18 minutes. Always ensure mechanical integrity before each use.

Home care use:

- Never use excessive force to insert or remove the dilator from the puncture. If you are unable to insert or remove the dilator with gentle force, immediately consult your clinician.

PRECAUTIONS

General:

- It is advisable to lubricate the dilator before the procedure. Only use water-soluble lubricants in an appropriate amount. Do not use any kind

of oil or ointment during the procedure. Oil based lubricants (e.g. Vaseline) will weaken, damage or destroy the device and shall not be used at any time.

- Make sure that the diameter of the dilator is adequate for the puncture. Over-dilatation may result in leakage around the voice prosthesis.
- Never use excessive force to insert or remove the dilator from the puncture, this may cause bleeding or a false route.
- Do not use a contaminated dilator. A contaminated dilator may cause infection and inflammation.
- Do not pull the safety medallion during removal. See Operating instruction.
- Do not use the dilator for a too long period of time as this may damage the esophagus (food pipe) wall.

Professional use:

- 
- Patients on anticoagulant therapy should be carefully evaluated for the risk of hemorrhage prior to dilatation of the TE puncture.
- 

Home care use:

To reduce the risk of infection:

- Make sure that your hands are thoroughly clean and dry before using your hands in the area of your tracheostoma, and before handling your voice prosthesis and/or accessories
- Clean and disinfect the device after each use according to instructions in section Cleaning and sterilization below.
- If the dilator looks dirty or has air dried in an area with a risk of contamination the dilator should be both cleaned and disinfected before use. A risk of contamination could exist if the device has been in contact with a pet, dropped on the floor or other

gross contamination, or if someone has a respiratory infection.

- During hospitalization it is important to both clean and disinfect the dilator after use but also before use since there is an increased risk of infection and contamination. In a hospital it is better to use sterile water for cleaning and rinsing, rather than drinking water.

Adverse events

Consult your clinician if you:

- bleed from the puncture tract. Bleeding may be caused by damage of the puncture tract, damage of the esophageal wall or by bleeding from esophageal vessels. Minor bleeding during dilatation is normal. If a minor bleeding does not stop, you should consult your clinician.
 - accidentally dislodge the dilator. The dilator may be used for stenting the puncture during a longer period of time. Accidental dislodgement during these circumstances may or may not be followed by shrinkage of the puncture and subsequent closure.
 - experience an allergic reaction against the lubricating (anaesthetizing) gel. Read the instructions for use accompanying the used lubricant.
- 
- 

Instructions for use

Checkout procedures

General:

Always handle the dilator with clean hands or gloves and avoid touching parts that will enter the puncture.

Prior to each use check the integrity of the safety strap by slightly pulling the medallion with one hand while holding the device handle with the other hand (Figure 2). Discard the device

if it shows any sign of damage.

It is advisable to lubricate the dilator with some water soluble (anaesthetizing) gel prior to use.

Operating instruction

Home care use:

Make sure that you have a magnifying mirror and sufficient light to illuminate the tracheostoma.

Avoid swallowing once the voice prosthesis is removed and the puncture is open since this may lead to aspiration of saliva and/ or coughing. Keeping your mouth open can help you avoid swallowing for short periods.

Do not clean, remove, or otherwise manipulate the Provox NID voice prosthesis (or use the Provox Dilator) immediately after a meal since the stimulation of the puncture tract immediately after eating may lead to gagging and/or vomiting.

Make sure that you have completely understood and fully practiced all procedures under the supervision of your clinician before you perform any procedure on your own.

Rinse the dilator in drinking water before use (Figure 3).

Lubricate the part of the dilator that will enter the puncture using a thin layer of water soluble (anaesthetizing) gel (Figure 4) if required.

Remove the Provox NID voice prosthesis from the TE puncture before inserting the dilator. Make sure the puncture tract is empty and/or open.

Probe the puncture with the tapered tip of the dilator with the curve of the dilator directed downwards, and slowly and gently insert the device until you reach the stop in front of the handle (Figures 5-7).

Keep the dilator in this position for the

desired amount of time. The time needed to sufficiently dilate the puncture tract to facilitate voice prosthesis insertion varies from person to person.

You shall carefully tape the safety medallion to the skin while you wait (Figure 8).

Once the dilator has been in the puncture tract for the desired amount of time, gently remove the dilator by pulling the handle (Figure 9). Do not pull the safety medallion.

Insert the Provox NID voice prosthesis using the standard procedure.

Problem solving:

There might be several possible causes to why you are unable to insert the dilator, such as wrong angle/direction when inserting, shrinkage of puncture or deformation of the puncture tract. If this should occur never use force to insert the dilator, you must consult your clinician immediately.

If you have inserted the dilator too far into the puncture tract, you might have used the product in the wrong way or might have an enlarged puncture tract. If this occurs gently try to retrieve the device by pulling the medallion. You should never use force. If you cannot retrieve the device, consult your clinician immediately.

Professional use:

Remove the voice prosthesis from the TE puncture.

Probe the puncture with the tapered tip of the dilator of adequate diameter, with its curvature directed downwards. Gently insert the device until the stop of the handle.

Keep the dilator in this position for the desired amount of time. The waiting time needed for dilation can vary. Tape the medallion to the skin during the waiting time.

Cleaning and sterilization

Home care use:

Clean and disinfect the dilator after each use.

To clean the dilator mix 2 drops of dish soap in 250 ml (8,5 fl oz) warm, 50°-100°C/122°-212°F, drinking water. Place the dilator in the mixture (Figure 10).

Gently rub the dilator for 15 seconds to remove mucus and debris (Figure 11).

Rinse it in drinking water (Figure 12).

Disinfect the dilator by soaking it in hydrogen peroxide 3% for 60 minutes or ethanol 70% for 10 minutes or isopropyl alcohol 70% for 10 minutes (Figure 13).

After cleaning and disinfecting the device, place the dilator on clean gauze and leave it to air dry (Figure 14). When the dilator is dry, store it at room temperature in a clean closed container. Make sure that the stored device is protected from both heat and direct sunlight.

Caution: Do not clean or disinfect by any other method as this might cause product damage and injury.

Professional use:

Cleaning, disinfection and sterilization are described in “Instructions for cleaning and sterilization”. If this information is missing, please contact your supplier or visit our website on www.atosmedical.com.

The material in the Provox Dilator withstands 50 cleaning and steam sterilization procedures at 137°C/278°F for 3 minutes or at 134°C/273°F for 18 minutes. Always ensure mechanical integrity before each use.

How supplied

The device is manufactured in a cleanroom and supplied non-sterile.

Device lifetime

Caution: The dilator(s) shall be replaced every 2 years, or after 50 sterilization cycles, whichever comes first. Check mechanical integrity of the device prior to each use according to the Checkout procedures.

Disposal

After use, the product may be a potential biohazard. Handling and disposing should be carried out in agreement with medical practice and applicable national laws and legislations.

Ordering information

See last page of this manual.

Date of printing

See version number on the back cover of this manual.

User assistance information

For additional help or information please see back cover of this manual for contact information.

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© Atos Medical AB, Sweden. 201201A, REF 10322