
# PROVOX<sup>®</sup> Micron HME<sup>™</sup>

## Instructions for Use

**Rx  
ONLY**

**CE**

# Atos

Atos Medical **Your voice**

Figure 1


Figure 2


Figure 3



Figure 4

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and copyrights**

Provox<sup>®</sup> is a registered trademark owned by Atos Medical AB, Sweden.

Provox<sup>®</sup> Micron HME<sup>™</sup> is a trademark of Atos Medical AB. For information about patents, please refer to the webpage [www.atosmedical.com](http://www.atosmedical.com).

# Contents

ENGLISH .....	5
DEUTSCH .....	8
NEDERLANDS.....	11
FRANÇAIS .....	14
ITALIANO.....	16
ESPAÑOL .....	19
PORTUGUÊS.....	22
SVENSKA .....	25
DANSK.....	27
NORSK.....	30
SUOMI .....	33
MAGYAR .....	35
POLSKI.....	38
ΕΛΛΗΝΙΚΑ .....	41
РУССКИЙ .....	44
日本語 .....	47

The Instructions for Use, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

## Intended use

The Provox<sup>®</sup> Micron HME<sup>™</sup> is a heat and moisture exchanger (HME) and air filtration device for laryngectomized patients. Provox Micron HME partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Provox Micron HME is intended to be used with the attachment devices in the Provox HME System.

## CONTRAINDICATIONS

This device shall not be used by patients who are unable to handle or remove the device themselves when needed, unless the patient is under constant supervision of a clinician or a trained caregiver. For example: patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at a risk for unpredictable periodic loss of consciousness.

## Description of the device

Provox Micron HME is a Heat and Moisture Exchanger combined with an electrostatic filter. The HME is a foam which contains a salt (Calcium Chloride). The HME retains the heat and moisture of the exhaled air. When inhaling, the retained heat and moisture in the HME is given back to the lungs. HME use may help improve the function of the lungs and reduce problems with e.g. coughing and mucus production. New users may experience slight discomfort in the beginning, related to increased breathing resistance. During the first weeks of use, mucus production may seem to increase. This is normal and means that the mucus is getting thinner and easier to cough up.

After a few weeks of HME use, this should stabilize and coughing and mucus production usually decreases.

The Provox Micron HME lid can be pressed down to occlude the stoma in order to speak with a voice prosthesis. When the pressure is released, the lid automatically comes up and the airway passage opens.

Provox Micron HME helps to filter inhaled air through consistent normal use. Thereby, small airborne particles, e.g., bacteria, viruses, dust and pollen are restricted from passing through the device into the lungs (see technical data below).

**Note:** Provox Micron is not intended to be used as a Personal Protective Equipment during work that requires breathing protection.

## Technical data

Pressure drop after 1 h at: (According to ISO 9360)	30 l/min	0.7 hPa
	60 l/min	1.9 hPa
	90 l/min	3.5 hPa
Pressure drop after 24 h at: (According to ISO 9360)	30 l/min	0.7 hPa
	60 l/min	1.8 hPa
	90 l/min	3.5 hPa
Moisture output: (According to ISO 9360)	22 mg/L	
Moisture loss at VT=1000 ml: (According to ISO 9360)	22 mg/L	
Bacterial Filtration Efficiency (BFE): (According to ASTM F2101)	≥ 99 %	
Virus Filtration Efficiency (VFE): (According to MIL-M-36954C and ASTM F2101)	≥ 99 %	

## WARNING

Unintentional or accidental pressure on the Provox Micron HME lid may cause difficulty in breathing or suffocation.

Inform the patient, caregivers and others of this feature to ensure that they understand the closing function. By putting pressure on the Provox Micron HME lid, the airway passage is closed. Closing the airway in order to allow voicing is a well-known feature for the laryngectomized patient with a voice prosthesis, but may be unknown for patients without a voice prosthesis.

## PRECAUTIONS

Provox Micron HME provides good protection through consistent normal use, as long as there is no air leakage. However, since there are other pathways for e.g. viruses and bacteria to enter the human body, total protection can never be guaranteed.

The same device must not be used for more than 24 hours after initial use. This can increase the risk for infection due to growth of e.g. bacteria.

Do not wash and re-use the device. Washing the HME impairs the filtering- and HME functions.

Do not disassemble Provox Micron HME. Disassembly will destroy its function.

Do not administer medicated nebulizer treatment over the device since the medication can be deposited in the device.

Do not use humidifiers or heated humidified oxygen over the device since the HME will become too wet.

## Instructions for use

Provox Micron HME can easily be attached (Fig. 1, Fig. 2) and removed (Fig. 4) from the Provox HME System attachment devices when needed. In order to get good protection, make sure that the seal is airtight by closing the Provox Micron HME and check for leakage.

In order to speak with a voice prosthesis, press down the lid (Fig. 3). This will direct the exhaled air through the voice prosthesis. When the pressure is released the airflow will pass through Provox Micron HME again.

Provox Micron HME can be removed during or after coughing if the stoma needs to be cleaned from mucus.

If you experience any problems with the product, contact your clinic.

## Device lifetime

Replace Provox Micron HME when needed. To ensure its proper function, the same device must not be used for more than **24 hours** after initial use. Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## Date of printing, Ordering and User assistance information

The version number of this Instruction for Use and manufacturer, contact information is located on the back cover. For ordering information, see end of this Instructions for Use.

## Attachment devices

The Provox Micron HME and other Provox HME Cassettes are intended to be used with attachment devices and accessories in Provox HME System. For more information read Instructions for Use for each of the products.

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© Atos Medical AB, 2018.