
# PROVOX<sup>®</sup>

## GuideWire

**Rx**  
ONLY

**CE** 0413

**STERILE EO**

# Atos

Atos Medical **Your voice**

# CONTENTS

ENGLISH.....	8
DEUTSCH .....	13
NEDERLANDS.....	18
FRANÇAIS .....	23
ITALIANO.....	28
ESPAÑOL .....	33
PORTUGUÊS.....	38
SVENSKA .....	43
DANSK.....	48
NORSK.....	53
SUOMI .....	58
ČESKY .....	63
POLSKI.....	68
HRVATSKI .....	73
TÜRKÇE .....	78
РУССКИЙ .....	83
英語.....	89
한국어.....	94
繁體中文 .....	99
עברית.....	107
العربية.....	112
ORDERING INFORMATION .....	113

Fig.1


Fig. 2







Fig. 3










Fig. 4

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and trademarks**

Provox® is a trademark owned by Atos Medical AB, Sweden.

The Instructions for Use, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

## 1. Descriptive Information

### 1.1 Indications for use

The **Provox GuideWire** is a sterile single use insertion device intended for placement of a sterile Provox indwelling Voice Prosthesis after total laryngectomy (primary or secondary puncture), or for retrograde replacement of a sterile Provox indwelling Voice Prosthesis.

### 1.2 Contraindications

Do not use if the patient has anatomical abnormalities, e. g. significant pharyngeal stenosis above the puncture site or severe trismus. Significant pharyngeal stenosis may preclude insertion of the voice prosthesis. Severe trismus may preclude proper protection of the pharyngeal wall during secondary puncture leading to harm of the esophageal tissue.

### 1.3 Description of the device

The Provox GuideWire is a sterile, single use device for introduction and replacement of sterile Provox indwelling Voice Prostheses. The GuideWire has a connector for attachment of the safety strap of the new voice prosthesis and an 8 mm Stopper for transoral removal of the remnant of the old voice prosthesis.

The Provox GuideWire package contains the following items:

- 1 Provox GuideWire, sterile (Fig.1),
- 1 Instructions for use Provox GuideWire, non-sterile.

The Provox GuideWire has two functional parts (Fig.1)

1. Stopper – for transoral removal of the old Provox Voice Prosthesis (made of ABS plastic) (Fig.1.1),
2. Connector – for attachment of a Provox Voice Prosthesis safety strap (made of nylon plastic material (Polyamide) (Fig.1.2).

The Provox GuideWire tube is made of a PVC plastic material.

## 1.4 WARNINGS

### Pre-surgery

- **DO NOT** use the product if the package or product is damaged or opened. Unsterile product may cause infection.
- **DO NOT REUSE** and **DO NOT RESTERILIZE** by any method. Intended for single use only. Reuse may cause cross contamination. Cleaning and reesterilization may cause structural damage.

### During insertion

- **DO** always proceed slowly and without using excessive force when pulling the GuideWire through. Tissue damage or dislodgement or extrusion of parts of the Provox GuideWire may otherwise occur which can lead to subsequent ingestion or tissue damage.
- **DO NOT** proceed if the voice prosthesis gets stuck in pharynx or esophagus, this may lead to tissue damage.

For further information about these events, see section Adverse events.

## 1.5 PRECAUTIONS

Always assess the suitability before using the Provox GuideWire in the patient. In cases with lacking suitability, do not use the product.

- **DO** carefully assess patients with bleeding disorders or patients undergoing anticoagulant treatment for the risk of bleeding or hemorrhage prior to placement or replacement of the voice prosthesis.
- **DO** ensure that in case of known or suspected pharyngeal stenosis, dilatation of the pharynx should be performed before using the GuideWire.
- **DO** always use aseptic technique during the procedure in order to reduce the risk of infection.
- **DO** ensure that the safety strap of the voice prosthesis is adequately attached to the GuideWire. If proper attachment is not achieved, the voice prosthesis may come loose from the GuideWire, causing failure to complete the procedure.
- **DO NOT** load or reload the voice prosthesis if the safety strap has been cut off or is broken.

## 1.6 Adverse events and troubleshooting information

**Aspiration of parts of the Provox GuideWire or parts of the Provox System** – Accidental aspiration of parts of the Provox GuideWire or parts of the Provox System may occur. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Immediate symptoms may include coughing, wheezing or other abnormal breathing sounds, dyspnea, and respiratory arrest, partial or inadequate air exchange and/or asymmetrical chest movement with respiration. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma. If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If aspiration of the object is suspected, the object should be located and retrieved endoscopically using a non-toothed grasping forceps.

**Ingestion of parts of the Provox GuideWire or parts of the Provox System** – Accidental ingestion of parts of the Provox GuideWire or parts of the Provox System may occur. As with any other foreign body, the symptoms caused by ingestion of parts of the GuideWire or parts of the Provox System, depend largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach. Foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies from the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs or the object fails to pass through the intestinal tract.

Spontaneous passage of the object may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested object. If the object does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted. The device may be retrieved by using a non toothed grasping forceps.

**The Provox GuideWire gets stuck in the pharyngeal mucosal wall** – The Provox GuideWire may get stuck in or be interfered by the pharyngeal mucosal wall. With slight pressure, the Provox GuideWire will generally bend near the tip and slide upwards towards the pharynx. Stop the procedure if the general bending does not help.

**Hemorrhage/Bleeding of the puncture** – Slight bleeding from the edges of the tracheoesophageal (TE) puncture may occur during intended use of Provox GuideWire and generally resolves

spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage prior to secondary puncture and prosthesis placement.

**Rupture of the tracheoesophageal tissue** – In case of a rupture of the TE tissue the TE puncture procedure should be abandoned and the rupture should be sutured immediately. TE puncture should only be repeated after proper healing of the tissues.

**The Provox GuideWire gets stuck inside the esophagus** – If the Provox GuideWire gets stuck inside the esophagus, do not use excessive force. Attempts to pull out a stuck GuideWire may cause tissue damage and/or breaking the GuideWire. The device has to be retrieved with an endoscope.

## 2. Instructions for use

### 2.1 Preparation

Check the integrity of the sterile package. Do not use the product if the package is damaged or opened.

### 2.2 Operating instruction

#### 2.2.1 Primary puncture placement (Fig 2):

1. Insert the Provox Pharynx Protector in the open pharynx (Fig 2.1).
2. Use the Provox Trocar with Cannula to create a TE-puncture (about 8-10 mm from the edge of the tracheostoma) (Fig 2.2).
3. Remove the Trocar after the puncture, while leaving the Trocar Cannula in place (Fig 2.3).
4. Introduce the Provox GuideWire in anteroposterior direction through the opening of the Trocar Cannula.

The connector of the GuideWire appears out of the upper opening of the Pharynx Protector, which then can be removed (Fig 2.4).

**CAUTION:** Do not attach the voice prosthesis when the Pharynx Protector is in place, since the voice prosthesis can get stuck.

5. The safety strap of the voice prosthesis of the proper model, diameter and size is attached and secured to the GuideWire connector. The Trocar Cannula is removed from the puncture opening and is pulled back.

Careful traction on the GuideWire will bring the strap of the prosthesis into the TE-puncture. The voice prosthesis is carefully rotated into position with two non-toothed hemostats firmly holding the tracheal flange (Fig 2.5).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

6. When the voice prosthesis is correctly positioned, the safety strap is cut off while still attached to the GuideWire and the voice prosthesis is turned so that the part where the safety strap was attached points downwards into the trachea. The laryngectomy procedure can then be completed (Fig 2.6).

#### 2.2.2 Secondary puncture placement (Fig 3):

1. Insert a rigid endoscope (or equivalent instrument) into the esophagus (Fig 3.1).
2. Use the Provox Trocar with Cannula to create a TE-puncture (about 8-10 mm from the edge of the tracheostoma) (Fig 3.2).
3. Remove the Trocar after the puncture, while leaving the Trocar Cannula in place (Fig 3.3).
4. Introduce the Provox GuideWire through the Trocar Cannula. The GuideWire is pushed upwards through the scope until it is stopped by the Stopper at the end of the Cannula. The endoscope can then be removed (Fig 3.4).

**CAUTION:** Do not attach the voice prosthesis when the endoscope is in place, since the voice prosthesis can get stuck.

5. The safety strap of the voice prosthesis of the proper model, diameter and size is attached and secured to the GuideWire connector. The Trocar Cannula is removed from the puncture opening and is pulled back.

Careful traction on the GuideWire brings the strap of the voice prosthesis into the TE-puncture. The voice prosthesis is carefully rotated into position with two non-toothed hemostats firmly holding the tracheal flange (Fig 3.5).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

6. When the voice prosthesis is correctly positioned, the safety strap is cut off while still attached to the GuideWire and the voice prosthesis is turned so that the part where the safety strap was attached points downwards into the trachea (Fig 3.6).

### **2.2.3 Replacement procedure (Fig 4):**

#### **Anterograde vs. retrograde insertion**

If the sterile indwelling Voice Prosthesis cannot be inserted in an anterograde manner, retrograde insertion by use of the Provox GuideWire may be performed.

Replacement should be carried out as an outpatient procedure using local anaesthesia. Occasionally, the procedure may be carried out under general anesthesia upon determination by the physician.

#### **Anterograde removal of the old prosthesis (not Provox(1)):**

1. Grasp the tracheal flange with a non-toothed hemostat and pull the prosthesis out from the TE puncture.
2. Insert the GuideWire through the TE puncture and push the GuideWire upwards through the pharynx and out of the mouth.

#### **Retrograde removal of the old prosthesis (Only Provox(1)):**

The old Provox(1) voice prosthesis can be removed transorally by the use of the Provox GuideWire.

1. While the old voice prosthesis is still in situ, introduce the GuideWire through the old prosthesis and push the GuideWire upwards through the pharynx and out of the mouth.
2. Grasp the tracheal flange of the prosthesis with a non-toothed hemostat. Cut it off the tracheal flange with a scalpel and remove it. With the tracheal flange removed, the Stopper of the GuideWire will pull the remains of the old voice prosthesis into the esophagus and up through the mouth. Discard the old voice prosthesis.

#### **Placement of a sterile indwelling Provox Voice Prosthesis:**

The GuideWire Connector is now outside the mouth (Fig 4.1).

3. Attach and secure the new voice prosthesis to the GuideWire Connector. Pull the GuideWire downwards, towards the existing TE-puncture. Careful traction on the GuideWire brings the strap of the voice prosthesis in to the TE-puncture (Fig 4.2).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

4. Introduce the tracheal flange of the voice prosthesis into the puncture with the help of a non-toothed hemostat. Then turn the voice prosthesis so that the safety strap is pointing downwards, into the trachea. Let the patient drink and swallow while observing the voice prosthesis, to ensure no leakage occurs. Cut off the safety strap while it is still attached to the GuideWire and discard the GuideWire (Fig 4.3).

## **2.3 Cleaning and sterilization**

The Provox GuideWire is provided sterile (EO) and intended for single use only and can NOT be cleaned or resterilized.

## **2.4 Disposal**

After use, the product may be a potential biohazard. Handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations.

## 2.5 Additional devices

**Provox Trocar with Cannula:** Stainless steel instruments for performing a primary or secondary TE puncture. The Provox Trocar creates a small, rounded puncture, while the Cannula provides a guided passage for the GuideWire through the freshly made TE puncture.

**Provox Pharynx Protector:** Stainless steel instrument for protection of the pharyngeal wall during primary TE puncture.

## 3. Additional information

### 3.1 Ordering information

See back of this Instruction for Use.

### 3.2 Date of printing

See version number on the back cover of this Instruction for use.

### 3.3 User assistance information

For additional help or information please see back cover of this manual for contact information.

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© Atos Medical AB, 2016.