

PROVOX[®] Vega[™]

Puncture Set

MD

Rx ONLY

CE 0413

STERILE EO

42° C
108° F
RT
2° C
36° F

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden. Provox® Vega™ is a trademark of Atos Medical AB. For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Producent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamin-tojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιητής; Ístehsalçı; Производитель; Produsen; Pengilang; 製造業者; 제조사; 製造商; 製造商; היצרן; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Ἰσθημ τήρησι; მწარმოების თარიღი; Ἀρτοποιηման ամսաթիվ; Ístehsal tarixi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조일자; 製造日期; 製造日期; תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Användes inden; Brukes innen; Viimeinen käyttöpäivä; Síðasti notkunardagur; Kölblík kuni; Derīguma termiņš; Naudoti iki; Datum expirace; Felhasználhatóság dátuma; Dátum najneskoršej spotreby; Rok uporabnosti; Data ważności; Data expirării; Rok uporabe; Rok upotrebe; Ημερομηνία λήξης; Срок на годност; Son kullanim tarihi; ვარგისობა; Ἰσθημ անվտանգობის ժամկետ; Son istifadə tarixi; Срок годности; Gunakan sebelum tanggal; Tarikh guna sebelum; 使用期限; 사용 기한; 使用截止日期; 使用截止日期; תאריך אחרון לשימוש; تاريخ انتهاء الصلاحية

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Please see accompanying Illustration manual for illustrations referenced in this Instructions for Use.

Stomal airway management

Emergency situations

It is important that emergency personnel knows that you are a neck breather. This will help them determining that oxygen and rescue breathing need to be administered to your stoma and not to your mouth or nose.

We recommend that you and your clinician download the guidelines for rescue breathing from www.atosmedical.com/rescuebreathing

Elective situations

If you need to undergo a procedure that requires intubation (putting a breathing tube in your windpipe), it is very important that the anesthesiologist and the doctor who is performing the procedure are aware that you are a neck breather and that you are using a voice prosthesis. It is important for them to understand that the voice prosthesis should stay in place. If it is removed, fluids from your mouth, food pipe, or stomach may get into your windpipe.

It is very important that the intubation tube is inserted and removed carefully so that the voice prosthesis does not dislodge or come out.

1. Descriptive information

1.1 Indications for use

Provox Vega Puncture Set is a device for performing a primary or secondary tracheo-esophageal (TE) puncture in laryngectomized patients, with integrated placement of a Provox Vega voice prosthesis.

The Provox Vega voice prosthesis is a sterile single use indwelling voice prosthesis intended for voice rehabilitation after surgical removal of the larynx (laryngectomy). Cleaning of the voice prosthesis is performed by the patient while it remains in situ.

1.2 CONTRAINDICATIONS

Do not use the Provox Vega Puncture Set if the patient has anatomical abnormalities that may hinder safe puncturing of the TE wall or safe voice prosthesis placement (e.g., significant stenosis or significant fibrosis at the puncture site) as this may cause tissue damage.

Do not use the Provox Vega Puncture Set for secondary TE puncture if the patient suffers from severe trismus that precludes proper protection of the pharyngeal wall. Failure to protect the pharynx during puncture may lead to unintended trauma of the pharyngeal/ esophageal tissue.

1.3 Description of the device

The Provox Vega Puncture Set is a device for creating a primary or secondary TE puncture, with subsequent dilatation of that puncture to a width that facilitates placement of the included Provox Vega voice prosthesis. The Provox Vega voice prosthesis is preloaded on the Puncture Dilator, which is part of the device.

The Provox Vega Puncture Set is intended for single use only and the package contains the following sterile items in a blister package (Fig. 1):

- 1 Pharynx Protector (Fig. 1.1) made of transparent thermoplastic,
- 1 Puncture Needle (Fig. 1.2) made of surgical stainless steel,
- 1 Guidewire (Fig. 1.3) made of pre-colored fluoroplastic,
- 1 Puncture Dilator with 1 preloaded Provox Vega voice prosthesis (Fig. 1.4). The Puncture Dilator is made of thermoplastic elastomer and polypropylene; and the Vega voice prosthesis is made of medical grade silicone rubber and fluoroplastic.

The preloaded Puncture Dilator contains the following functional characteristics and components:

- a dilator (Fig. 1.4.1),
- a dilator strap (Fig. 1.4.2) connecting the dilator and the voice prosthesis interface,
- a dilator loop (Fig. 1.4.3) constituting the voice prosthesis interface,
- a Wirelock (Fig. 1.4.4) containing interfaces to the voice prosthesis safety strap (Fig. 1.4.6) and to the Guidewire,
- a Provox Vega voice prosthesis (Fig. 1.4.5) with its safety strap (Fig. 1.4.6) connected to the Wirelock and oriented so that the tracheal flange (Fig. 1.4.7) of the voice prosthesis is facing towards the Wirelock.

The Provox Vega voice prosthesis contains a one-way valve that keeps the TE-puncture open for speech, while reducing the risk of fluids and food entering the trachea.

The Provox Vega voice prosthesis is not a permanent implant, and needs periodic replacement. The prosthesis is available in different diameters and several sizes.

The set also includes the following non-sterile items:

- 1 Instructions for use - Provox Vega Puncture Set (including 1 Illustration manual)
- 1 Provox Vega Patient's Manual
- 1 Provox Brush of a size corresponding to the voice prosthesis
- 1 Provox Brush Instructions for Use

1.4 WARNINGS

Pre-surgery

- **DO NOT** use the product if the package is damaged or opened. Unsterile product may cause infection.
- **DO NOT REUSE** and **DO NOT RESTERILIZE** by any method. This device is intended for single use only. Reuse may cause cross contamination. Cleaning and resterilization may cause structural damage to the device.
- **DO** proceed with great care if the patient has received radiotherapy with or without concurrent chemotherapy. These circumstances increase the risk of puncture-related complications (e.g., widening, granulation, atrophy). Therefore, ensure that the tissue integrity is sufficient for creation of a TE puncture.

During surgery

General

- **DO** ensure that the Pharynx Protector is inserted deep enough into the esophagus before performing primary TE puncture, by palpating the TE wall. Puncture without proper positioning of the Pharynx Protector may cause tissue damage.
- **DO** ensure that the Guidewire is inserted through the needle and through the lumen of the Pharynx Protector so that it does not damage the TE wall.
- **DO** ensure that the Puncture Needle is removed before removing the pharynx protector. The Puncture Needle may cause tissue damage without proper pharynx protection.
- **DO** ensure that the Puncture Dilator is mounted on the esophageal end of the Guidewire and not the tracheal. The dilatation should be performed in the posteroanterior direction. Dilatation in the wrong direction causes reversed placement of the voice prosthesis, which causes aspiration and inability to speak.
- **DO NOT** withdraw the Guidewire back through the Puncture Needle. Damage, shearing and/or scraping of the Guidewire may occur. If the Guidewire must be withdrawn, remove the Guidewire and the Puncture Needle simultaneously, as one unit, to prevent the Puncture Needle from damaging the Guidewire.
- **DO NOT** use toothed hemostats or other instruments that may damage the product.

Secondary puncture

- **DO NOT** use the included Pharynx Protector during Secondary punctures. It is only intended to be used during Primary punctures.
- **DO** ensure that the pharyngeal/esophageal tissue is adequately protected, e.g., by use of a rigid endoscope before performing secondary TE puncture.

Post surgery

Use of the voice prosthesis

Dislodgement or extrusion of the Provox Vega voice prosthesis from the TE puncture and subsequent ingestion, aspiration or tissue damage may occur. For further information about these events and how to prevent them, see section Adverse Events and Trouble Shooting Information below.

To reduce the risk of dislodgment or extrusion and its potential consequences:

- **DO** select the proper prosthesis size (i.e. length). A tight fit due to a too short voice prosthesis may cause tissue necrosis and extrusion.
- **DO** instruct the patient to only use genuine Provox accessories of corresponding size and diameter (e.g., Brush, Flush, Plug) for maintenance and to avoid all other kinds of manipulation.
- **DO** instruct the patient to consult a physician immediately if there are any signs of tissue edema and/or inflammation/infection.
- **DO** choose laryngectomy tubes or stoma buttons, if used, with a suitable shape that do not exert pressure on the prosthesis during use, or catch onto the tracheal flange of the prosthesis during insertion and removal of the tube or button.

1.5 PRECAUTIONS

Always assess the suitability of the tissue in the area of the TE puncture. In cases with lacking suitability, e.g. due to excessive scar tissue or radiation fibrosis, proceed with great care and abort the procedure if dilatation of the TE puncture requires too much force.

- **DO** carefully assess patients with bleeding disorders or patients undergoing anticoagulant treatment for the risk of bleeding or hemorrhage, prior to secondary puncture and prosthesis placement.
- **DO** always use aseptic technique when handling the Puncture Set in order to reduce risk for infection.
- **DO** remove the Pharynx Protector before initiating dilatation. The voice prosthesis may get stuck inside the Pharynx Protector if attempting to complete the procedure without removing the Pharynx Protector.
- **DO** ensure that the Guidewire is adequately threaded and locked in position in the Wirelock. If proper locking is not achieved, the Guidewire may come loose from the Wirelock, causing failure to complete procedure.
- **DO** always proceed slowly and without using excessive force during dilatation and placement of the prosthesis. Tissue damage may otherwise occur.
- **DO** support the TE tissue during dilatation. Otherwise, rupture may occur. In case of a rupture of the TE tissue the TE puncture procedure should be abandoned and the rupture should be sutured immediately. TE puncture should only be repeated after proper healing of the tissues.
- **DO NOT** reload after the safety strap of the prosthesis has been cut off, as the safety mechanism in that case has been compromised, with risk for dislocation of the voice prosthesis during the procedure.

1.6 Adverse events and troubleshooting information

1.6.1 During use of the Puncture Set

(Sub-) mucosal injury

During puncture, the Puncture Needle or the Guidewire could cause a (sub) mucosal injury if the Pharynx Protector is not located in its correct position, or if the puncture is made incorrectly. In case of suspicion of (sub)mucosal damage, it is recommended that the patient receives a prophylactic course of postoperative antibiotics and is fed through a nasogastric tube or similar to allow healing.

Re-puncture

If the Pharynx Protector is not positioned appropriately when puncturing, or if, for any other reasons, there is need for re-puncture, the tissue should be assessed for possible sub-mucosal injury and the puncture procedure repeated with the Pharynx Protector positioned correctly.

Forgot to remove the Pharynx Protector

If the Puncture Dilator is connected to the Guidewire prior to removal of the Pharynx Protector, the procedure cannot be completed. If this occurs, detach the Puncture Dilator from the Guidewire and remove the Pharynx Protector. Also see Reloading below.

Reloading

In some cases, reloading of the Puncture Dilator may be necessary; e.g., if the voice prosthesis is pulled completely through the puncture during prosthesis placement. The reloading procedure is described in section 2.2.3 and illustrated in Fig. 4.

1.6.2 During use of the voice prosthesis

Dislodgement of the voice prosthesis– Dislodgement can be caused by infection and/or edema of the TE-puncture, granulation around the puncture or hypertrophic scarring around the puncture. Dislodgement might lead to aspiration or ingestion of the prosthesis. See below.

Aspiration of the prosthesis – Accidental aspiration of the voice prosthesis or other components of the voice rehabilitation system may occur. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Immediate symptoms may include coughing, wheezing or other abnormal breathing sounds, dyspnea, and respiratory arrest, partial or inadequate air exchange and/or asymmetrical chest movement with respiration. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma.

If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If aspiration of the device is suspected, a CT scan of the lungs should be performed to confirm aspiration and locate the device. If the CT scan confirms aspiration of the device, the device may be retrieved endoscopically using a non-toothed grasping forceps.

The silicone housing of the Vega voice prosthesis can also be located endoscopically. On a CT scan and during endoscopy, the device may appear as an oval shape with an opening in the middle with an outer diameter of about 10-17 mm (the flanges of the device), or as a cufflink shape with a shaft length of 8, 10, 12.5 or 15 mm, depending on the size of the device. During endoscopy, reflections from the light source on the clear silicone rubber may be seen. Also, in prostheses that have been in situ for some time, white or yellow appearing Candida deposits may be visible on the device.

Ingestion of the prosthesis – Accidental ingestion of the voice prosthesis, or other components of the voice rehabilitation system, may occur. As with any other foreign body, the symptoms caused by ingestion of the prosthesis or a component of the voice rehabilitation system depends largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the lower esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach. Foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies from the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs or the object fails to pass through the intestinal tract.

Spontaneous passage of the device may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested device. If the device does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted. The silicone housing of the Vega voice prosthesis can be located and retrieved endoscopically. The device may be retrieved by using a non toothed grasping forceps. During endoscopy, the device may appear as an oval shape with an opening in the middle with an outer diameter of about 10-17 mm (the flanges of the device), or as a cufflink shape with a shaft length of 8, 10, 12.5 or 15 mm, depending on the size of the device. Reflections from the light source on the clear silicone rubber may be seen. In prostheses that have been in situ for some time, white or yellow appearing Candida deposits may be visible on the device.

Infection and/or edema of the TE-puncture – Infection, granulation formation and/or edema of the puncture may increase the length of the puncture tract. This may cause the prosthesis to be drawn inward and under the tracheal or esophageal mucosa. Inflammation or overgrowth of the esophageal mucosa may also cause the prosthesis to protrude from the puncture. Temporary replacement of the prosthesis by a prosthesis with a longer shaft is then advisable. If standard medical treatment does not resolve the infection, the prosthesis should be removed. In some cases stenting the puncture with a catheter might be considered. If the puncture closes spontaneously re-puncture for insertion of a new prosthesis may be required.

Granulation around the puncture – Formation of granulation tissue around the TE-puncture has been reported at an incidence of about 5%. Electrical, chemical, or laser cauterization of the area of granulation may be considered.

Granulation/Hypertrophic scarring around the puncture – Bulging of the tracheal mucosa over the tracheal flange may occur if the prosthesis is relatively short. This excess tissue may be removed by using a laser. Alternatively, a prosthesis with a longer shaft can be used.

Protrusion/extrusion of the prosthesis – Protrusion of the prosthesis and subsequent spontaneous extrusion is sometimes observed during infection of the TE-puncture. Removal of the prosthesis is required to avoid dislodgement into the trachea. The puncture may close spontaneously secondary to the removal of the prosthesis. Re-puncture may be necessary for insertion of a new prosthesis.

Tissue damage – If the prosthesis is too short, too long, or is pushed frequently against the esophageal wall by a laryngectomy tube, stoma button, or the patients' finger, damage of the puncture, tracheal and/or esophageal tissues may occur. Inspect the conditions regularly to avoid severe damage.

Leakage through the valve – Leakage through the prosthesis may occur because:

- Candida overgrowth near the valve seat and valve flap may lead to incomplete closure of the valve flap, which causes leakage through the device. This is a normal event in prosthetic voice rehabilitation and an indication to replace the voice prosthesis.
- Stronger negative pressure in the PE-segment occurs during swallowing. To investigate this, transstomal inspection of the valve flap of the prosthesis should be performed during swallowing.

Leakage around the prosthesis – Transient leakage around the prosthesis may occur and may resolve spontaneously. The most common reason is that the prosthesis is too long, which is solved by inserting a shorter prosthesis. If the problem is not solved by inserting the correct length prosthesis, other factors that may affect tissue integrity in the area of the TE puncture (for example gastroesophageal reflux or thyroid function) should be considered, evaluated, and treated. In enlarged punctures with reduced retention strength other conventional treatment methods such as the injection of fillers (e.g., collagen) or temporary removal of the voice prosthesis, should be considered. If the leakage around the voice prosthesis is intractable, more conservative measures, such as surgical closure of the puncture may be necessary.

2. Instructions for use

Please see accompanying Illustration manual for illustrations referenced in this Instructions for Use.

For a visual overview of the different procedures you can find links to video animations under the section headings below.

Caution: The videos do not replace nor do they set forth the complete contents of the Instructions for Use and/or Prescriber Information, and are not a substitute for reviewing the entire contents of the Instructions for Use. The videos are only intended to further enhance the understanding of the procedure after review of the Instructions for Use.

2.1 Preparation

Prior to the puncture always determine what size and diameter of voice prosthesis to use. The appropriate size and diameter depends on the anatomy of the patient, local medical practice and preference of the surgeon.

2.2 Operating instruction

Check the integrity of the sterile package. Do not use the product if the package is damaged or opened. Unsterile product may cause infection.

2.2.1 Primary puncture and prosthesis placement

www.atosmedical.com/primary-puncture

1. After removal of the larynx and creation of the tracheostoma, before closure of the pharynx, insert the Pharynx Protector in the open pharynx/esophagus (Fig. 2.1).
2. Verify the correct location for the TE puncture by palpating the inside of the trachea at the desired puncture site. The oblique front opening of the Pharynx Protector (or the slit on the upper side, depending on surgical technique) should be felt during palpation (Fig. 2.2).
3. Insert the Puncture Needle at the correct puncture site (about 8-10 mm from the edge of the tracheostoma) until the tip of the needle reaches the inner lumen of the Pharynx Protector (Fig. 2.3).
If an endotracheal tube is in situ, this tube should be removed if it obstructs proper dilatation and integral placement of the voice prosthesis.
4. Insert the Guidewire into the hub of the Puncture Needle. Push the Guidewire through the needle until it extends approx. 20 cm out from the lumen of the Pharynx Protector (Fig. 2.4).
WARNING: Always verify that the Guidewire comes out through the lumen of the Pharynx Protector. Otherwise there is a risk for (sub) mucosal damage and the procedure needs to be restarted (see Adverse Events and Troubleshooting Information as well as Instructions for Reload of the Puncture Set).
5. Remove the Puncture Needle (Fig. 2.5).
CAUTION: Always remove the needle before removing the Pharynx Protector. There is a risk for damaging the esophageal tissue otherwise.
6. Remove the Pharynx Protector. Only the Guidewire should remain in situ before continuing (Fig. 2.6).
7. Insert the Guidewire extending from the esophageal side into the narrow end of the Puncture Dilator and push the Guidewire through the Puncture Dilator until it extends approx. 10 cm through the Puncture Dilator exit hole (Fig. 2.7).
8. Grab the tip of the Guidewire and insert it in the hole next to the exit hole (Fig. 2.8).
9. Tighten the Guidewire by pulling it from the narrow end of the Puncture Dilator and verify that it is secured to the Puncture Dilator (Fig. 2.9).
10. Using a continuous, smooth motion; dilate the puncture site by carefully pulling the Guidewire through the puncture. During dilatation, support the TE tissue (for example with two fingers) to reduce dilatation force. For better control, firmly grasp the Guidewire close to the Puncture Dilator (Fig. 2.10).
CAUTION: Dilatation and integral placement of the voice prosthesis should be carried out in the anterior/caudal direction with limited lateral movement in order to limit the force applied to the TE wall.
11. In the same continuous, smooth motion, carefully pull the Guidewire, Puncture Dilator and the Puncture Dilator loop through the puncture. The Puncture Dilator loop folds the tracheal flange of the voice prosthesis as the loop is pulled over the flange and through the puncture. The tracheal flange unfolds in the trachea when the loop releases it (Fig. 2.11).
Stop pulling immediately when the tracheal flange is released by the Puncture Dilator loop. If the tracheal flange does not unfold completely, it can be rotated in place using two non-toothed hemostats.
12. Grasp the tracheal flange of the voice prosthesis with a non-toothed hemostat, turn the prosthesis in the correct position, and cut the safety strap (Fig. 2.12).

2.2.2 Secondary puncture and prosthesis placement

www.atosmedical.com/secondary-puncture-and-prosthesis-placement

The Pharynx Protector (Fig. 1.1) included in the Provox Vega Puncture Set is not used during secondary puncture.

1. Choose an instrument (e.g., a rigid endoscope) which can function as;
 - a. a protector when the needle is penetrating the TE wall and
 - b. a guide for the correct location of the TE puncture and
 - c. a means to facilitate a safe passage of the Guidewire when it is passed through the pharynx and out of the mouth.Introduce the instrument into the esophagus (Fig. 3.1). If an endotracheal tube is in situ, this tube should be removed if it obstructs proper dilatation and integral placement of the voice prosthesis.
CAUTION: Always make sure that the instrument selected for pharynx protection does not contain a lumen for safe passage of the Guidewire. Otherwise, there is risk for (sub) mucosal damage.
2. Verify the correct location of the instrument by palpating the trachea at the intended puncture site. For additional visual and/or transilluminated guidance, a flexible endoscope could be used (Fig. 3.2).
3. Insert the Puncture Needle at the correct puncture site (about 8-10 mm from the edge of the tracheostoma) until the tip of the needle reaches the inside wall of the instrument (Fig. 3.3).
4. Insert the Guidewire into the hub of the Puncture Needle. Push the Guidewire into the needle, up through the lumen of instrument until it extends approx. 20 cm out through the distal end of the instrument (Fig. 3.4).
WARNING: Always verify that the Guidewire comes out through the lumen of the instrument chosen for pharynx protection. Otherwise there is a risk for (sub) mucosal damage and the procedure needs to be restarted (see Adverse Events and Troubleshooting Information as well as Instructions for Reload of the Puncture Set).
5. Remove the Puncture Needle (Fig. 3.5).
CAUTION: Always remove the needle before removing the instrument. There is a risk for damaging the esophageal tissue otherwise.
6. Remove the instrument used for pharynx protection. Only the Guidewire should remain in situ before continuing (Fig. 3.6).
7. From the cranial side, insert the Guidewire into the narrow end of the Puncture Dilator and push the Guidewire until it extends approx. 10 cm through the Puncture Dilator exit hole (Fig. 3.7).
8. Grab the tip of the Guidewire and insert it in the hole next to the exit hole of the Puncture Dilator (Fig. 3.8).
9. Tighten the Guidewire by pulling it from the narrow end of the Puncture Dilator and verify that it is secured to the Puncture Dilator (Fig. 3.9).
CAUTION: Ensure that the Guidewire is tightly secured in the Wirelock of the Puncture Dilator. If the Guidewire is not securely locked into the Wirelock, the Guidewire could separate from the Puncture Dilator and the Puncture Dilator may end up in the esophagus requiring retrieval using additional instruments (e.g., laryngeal forceps).
10. Using a continuous, smooth motion; dilate the puncture site by carefully pulling the Guidewire through the puncture site until the thick end of the Puncture Dilator has passed the puncture. During dilatation, support the TE tissue (for example with two fingers) to reduce dilatation force.

For better control, firmly grasp the Guidewire close to the Puncture Dilator (Fig. 3.10).

CAUTION: Dilatation and integral placement of the voice prosthesis should be carried out in the anterior/caudal direction with limited lateral movement in order to limit the force applied to the TE wall.

11. In the same continuous, smooth motion, carefully pull the Guidewire, Puncture Dilator and the Puncture Dilator loop through the puncture. The Puncture Dilator loop folds the tracheal flange of the voice prosthesis as the loop is pulled over the flange and through the puncture. The tracheal flange unfolds in the trachea when the loop releases it (Fig. 3.11).

Stop pulling immediately when the tracheal flange is released by the Puncture Dilator loop. If the tracheal flange does not unfold completely, it can be rotated in place using two non-toothed hemostats.

12. Grasp the tracheal flange with a non-toothed hemostat, turn the prosthesis in the correct position and cut the safety strap (Fig. 3.12).

2.2.3 Reload of the Puncture Set

www.atosmedical.com/reload-of-the-puncture-set

CAUTION: Reload of the Provox Vega Puncture Set must not be performed if the voice prosthesis safety strap has been cut or damaged during the first placement attempt.

If the procedure needs to be restarted the Provox Vega Puncture Set can be reloaded.

In order to reload the system:

1. Push the Guidewire from the narrow end of the Puncture Dilator until the Guidewire loosens from the Wirelock (Fig. 4.1-4.2).
2. Pull the Guidewire through the Puncture Dilator (Fig. 4.3).
3. Reload the Vega voice prosthesis in the Puncture Dilator ring (Fig. 4.4).

WARNING: The safety strap and tracheal flange of the voice prosthesis must be oriented towards the Puncture Dilator strap and the Wirelock when it is inserted into the Puncture Dilator loop (Fig. 4.4).

4. If needed, the Guidewire can be straightened to facilitate reinsertion.
5. Proceed according to “2.2 Operating Instructions”.

2.3 Cleaning and sterilization

The surgical components in the puncture set, as well as the voice prosthesis, are provided sterile (EO) and are intended for single use only and can NOT be cleaned or resterilized.

After placement, the voice prosthesis requires regular cleaning by the patient while the prosthesis remains in situ (see prosthesis maintenance below).

2.4 Important patient information

Inform the patient that:

- After a prosthesis placement traces of blood may be found in the sputum.
- Occasionally, mild leakage through or around the prosthesis may occur in the first weeks after insertion of a prosthesis. This often resolves spontaneously and does not require immediate replacement of the prosthesis.
- Speech training sessions with a speech language pathologist are advised in order to acquire optimal voice sound, fluent speech, and optimal intelligibility.

Ensure that the patient understands to contact their clinician if:

- Any changes in the appearance of the material of the prosthesis or in the way it fits in the puncture occur.
- Leakage occurs during eating and/or drinking, and cleaning the prosthesis does not help. Provox Vega Plug can be used for temporary preventing leakage during eating and drinking until the device can be replaced.
- Speaking becomes difficult (requires higher effort) and cleaning does not help.
- There are any signs of inflammation or tissue change at or near the puncture tract.
- Bleeding or overgrowth of tissue around the device.
- Persistent pain or discomfort in the region of the Vega voice prosthesis.
- Chronic cough, trouble breathing, or blood in the mucus. These can be signs of a serious health conditions that requires medical attention.

Prosthesis maintenance:

CAUTION: Only use genuine Provox accessories that are intended for use with Provox Vega when cleaning the prosthesis.

The patient should clean the prosthesis at least twice a day and after each meal with the Provox Brush by inserting the brush into the prosthesis and moving it gently back and forth with a twisting movement. After the brush is removed it should be wiped with a gauze. The procedure can be repeated as often as necessary. For details and how to clean the brush see Instructions for Use that come with the Provox Brush. In addition to using the Provox Brush, the patient may also clean the Provox Vega voice prosthesis with the Provox Flush. The Provox Flush can be used with drinking water or air. For details and how to clean the Provox Flush see Instructions for Use that come with the Provox Flush.

Compatibility with antifungal substances

In most cases treatment with antifungal agents should not be indicated but may be considered as a preventive measure if excessive candida overgrowth occurs on the prosthesis.

An unknown variety of chemical substances may influence the material properties of the device. Therefore, the introduction of antifungal medicines or drugs directly to or within close proximity of the voice prosthesis should be carefully assessed.

Laboratory testing shows no negative influence on the function of the Vega voice prosthesis and components when using the following antifungal medications: Nystatin, Fluconazole and Mikonazol.

Cleaning and disinfection of the accessories

The accessories should be cleaned after each use and disinfected at least once a day according to their Instructions for Use. During hospitalization there is an increased risk of device contamination and patient infection. Therefore, during hospitalization, it is important to clean and disinfect the accessories immediately after use and again just before use, and rinse using sterile water, rather than tap water.

For details and how to clean the Provox Accessories see Instructions for Use that accompanies the Provox Accessories.

2.5 Device lifetime of the Provox Vega voice prosthesis

Depending on individual biological circumstances the device life varies, and it is not possible to predict the integrity of the device over a longer period of time. The material of the device will be affected by e.g., bacteria and yeast, and the structural integrity of the device will eventually deteriorate.

Candida overgrowth of the prosthesis occurs in almost all patients. Radiotherapy, salivary content and dietary habits may influence how rapidly Candida can affect the silicone material and cause leakage through the prosthesis or other incompetence of the valve. Also see; Prosthesis maintenance under section 2.4 Important patient information.

The prosthesis is not a permanent implant, and requires periodic replacement. Depending on individual biological circumstances the device life varies, and it is not possible to predict the integrity of the device over a longer period of time. The prosthesis and especially the silicone material of the device will be affected by e.g., bacteria and Candida, and the structural integrity of the device will eventually deteriorate.

Indications for replacement of the Provox Vega voice prosthesis include leakage through the valve, blockage of the prosthesis, bacterial and candidal overgrowth leading to degradation of the materials and/or excessive pressure needed to obtain speech. Other reasons for earlier replacement could include medical indications such as problems with the puncture tract. Also see section 1.6 Adverse events and troubleshooting information / 1.6.2 During use of the voice prosthesis.

2.6 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

2.7 Accessories

Provox accessories are designed to be safe and effective for use with the Provox Vega voice prostheses. Do not use other devices since they may cause patient harm or product malfunction.

Accessories to maintain the device function (for patient use)

Provox Brush/Provox Brush XL: is used by the patient to clean the inside of the prosthesis.

Provox Flush: An additional cleaning device for patient use that allows flushing of the prosthesis.

Provox Plug: A tool for patient use that temporarily blocks the voice prosthesis.

Please see www.atosmedical.com or contact your local distributor for more info.

3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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