
PROVOX[®]

XtraFlange[™]

Rx
ONLY

CE 0413




STERILE | **EO**




Atos

Atos Medical **Your voice**

Illustrations

Contents

ENGLISH	6
DEUTSCH	11
NEDERLANDS.....	17
FRANÇAIS	23
ITALIANO.....	29
ESPAÑOL	35
PORTUGUÊS.....	41
SVENSKA	47
DANSK.....	52
NORSK.....	58
SUOMI	63
ČESKY	69
MAGYAR	74
POLSKI.....	79
HRVATSKI	85
ΕΛΛΗΝΙΚΑ	90
TÜRKÇE	96
РУССКИЙ	101
日本語	107
한국어	112
繁體中文	117
עברית	125
العربية	130
ORDERING INFORMATION	131

Clinician's manual and prescriber information

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden. Provox[®] XtraFlange[™] is a trademark of Atos Medical AB

The manual, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

Indications for use

Provox® XtraFlange™ is a silicone washer intended to reduce periprosthetic leakage that is detected on patients using indwelling Provox voice prostheses.

Placement is performed by a medical doctor or a trained medical professional in accordance with local or national guidelines.

Contraindications

The Provox XtraFlange shall NOT be used:

- in patients in whom the tracheoesophageal (TE) puncture is too wide to ensure adequate retention of the Provox voice prosthesis. A too wide puncture may increase the risk of dislodgement and aspiration of the device and/or the voice prosthesis. The Provox XtraFlange will NOT increase retention of the prosthesis!
- on Provox NID or other voice prostheses of any other brands. It may increase the risk of dislodgement and aspiration of the Provox XtraFlange and/or the voice prosthesis.

Device description

Provox XtraFlange is a white silicone washer that is intended to be placed between the tracheal flange of the prosthesis and the tracheal mucosa. It provides an extra seal against periprosthetic leakage through the adherence of the thin silicone sheet to the tracheal mucosa.

The device is supplied sterile and is intended for single use only.

The Provox XtraFlange is put in place with the help of two non-toothed hemostats while the prosthesis remains in situ (for appropriate Provox XtraFlange sizes, see order information). In order to prevent accidental aspiration during placement the Provox XtraFlange has a safety medallion. After proper placement, the safety medallion is cut off.

The Provox XtraFlange is intended to stay in place until the voice prosthesis is removed. The washer has to be discarded when the prosthesis is replaced.

WARNINGS

- Instruct the patient to seek immediate medical attention if the washer and/or prosthesis dislodges from the puncture. Immediate symptoms of aspiration may include coughing, wheezing or difficulty breathing. A foreign body in the airways may cause severe complications and has to be removed by a physician.
- A too tight fit may cause tissue necrosis and extrusion of the prosthesis. The washer is neither intended nor designed to make a too long prosthesis shorter. Instruct the patient to consult their clinician immediately if there are any signs of tissue edema and/or inflammation/infection. Treatment may require antibiotic therapy, temporary insertion of a prosthesis with a longer shaft, removal of the prosthesis and/or surgical closure of the puncture (with later possibly a repuncture).
- Provox XtraFlange is for single use; i.e. the washer has to be discarded when the prosthesis is replaced. Do not reuse the washer. Re-use and re-processing may cause crosscontamination and damage to the device, which could cause patient harm.

PRECAUTIONS

- Provox XtraFlange is not intended to be placed on the esophageal side of the voice prosthesis, and it is not intended to be glued in place.
- After placement check the Provox XtraFlange for damages. If damaged, replace it.
- Do not place more than one washer on the prosthesis.
- Do not cut off the safety medallion before placement since this may cause aspiration of the device.
- Instruct the patient to be careful when cleaning the voice prosthesis and stoma, and when inserting devices such as LaryTubes and LaryButtons into the stoma. Provox XtraFlange may be accidentally removed and aspirated. If this occurs the patient must seek immediate medical attention.

Instructions for use

Placement of the device

CAUTION: Always ensure that the tracheo-esophageal (TE) puncture has adequate retention of the Provox voice prosthesis by careful pulling at the tracheal flange with a hemostat. Place the Provox XtraFlange only if sufficient retention is confirmed. The Provox XtraFlange does not provide extra retention!

1. Have two non-toothed hemostats and scissors prepared (fig. 1).
2. Mount Provox XtraFlange on the first non-toothed hemostat (fig. 2).
3. Use this hemostat to grasp the tracheal (outer) flange of the voice prosthesis, preferably at the longest side of the flange (fig. 3).
Do not pull on the prosthesis. This hemostat must not be removed until the procedure is completed.
4. Grasp Provox XtraFlange with the second non-toothed hemostat. Place the washer by sliding it over the tracheal (outer) flange of the voice prosthesis in a circular motion (fig. 4).
5. Ensure that the washer is placed correctly between the tracheal flange of the prosthesis and the tracheal mucosa (rotate the washer) (fig. 5).
6. After Provox XtraFlange is placed over the tracheal flange and safely in position, cut off the safety medallion with a pair of scissors (fig. 6).

Replacement or removal

Remove and discard Provox XtraFlange when the voice prosthesis is changed. For detailed removal procedures of the voice prosthesis, please refer to the manual of the specific voice prosthesis.

If the voice prosthesis is to be removed anterograde

The voice prosthesis is pulled out of the TE puncture with a non-toothed hemostat, which also grasps the Provox XtraFlange. This way, both the Provox XtraFlange and the voice prosthesis are removed simultaneously. Both are to be discarded after removal.

If the voice prosthesis is to be removed retrograde

CAUTION: If the voice prosthesis is removed retrograde, the Provox XtraFlange must be removed separately prior to retrograde removal of the voice prosthesis to avoid the risk of aspiration of the Provox XtraFlange.

Grasp the Provox XtraFlange with a non-toothed hemostat. Make sure that the hemostat has a firm grip on the Provox XtraFlange before removal. If the Provox XtraFlange is aspirated it may cause partial airway obstruction and infection and must be removed immediately.

Placement of a new Provox XtraFlange after replacement of the voice prosthesis

To apply a new washer on the newly inserted voice prosthesis first observe whether it is necessary. If there is periprosthetic leakage, apply a new washer according to section “Placement of the device” above.

For instructions on how and when to change the voice prosthesis, please refer to the manual for the specific voice prosthesis.

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* United States Federal law restricts the sale, distribution and use by, or by order of a physician or a licensed practitioner of devices identified herein with the symbol "Rx". The prescription-free availability of these products outside the United States may vary from country to country.

**The availability of these products outside the United States may vary from country to country.

Atos

Atos Medical **Your voice**

Atos Medical AB

Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden

Tel: +46 (0)415 198 00 • Fax: +46 (0)415 198 98

www.atosmedical.com • info@atosmedical.com

© Atos Medical AB, 2017.