

**PROVOX**<sup>®</sup>  
**Life** Breathe better,  
whatever you do



Provox Life<sup>™</sup>  
Adhesives for greater  
comfort and fit.

# Provox Life™ Standard

Standard is an everyday adhesive that provides a strong seal around your stoma. A soft, flexible and low-profile adhesive that forms a strong seal and provides a secure fit.



## How does Provox Life™ Standard Adhesive differ from Provox FlexiDerm?

The adhesive material is the same. The difference is a slightly increased adhesive coating for Provox Life™ Standard Adhesive which makes it resist higher air pressures, which may help it to stick better and reduce blow-through.

### Key features:

- ✓ Strong (yet gentle) seal for normal skin
- ✓ Redesigned clover-shaped adhesive for easy application and better fit
- ✓ SecureFit™ for security, comfort and compatibility across all Provox Life™ HMEs



QR Code, watch our online video on how to apply Provox Life™ Standard in 4 simple steps

# Provox Life™ Stability

Stability Adhesive is a firm everyday adhesive that provides solidity, especially if a stoma is deep. It is an ideal adhesive for those who speak frequently or speak hands-free.

Provox® StabiliBase or  
Provox® XtraBase



Provox Life™ Stability



## Has the finger lift tab for Provox Life™ Stability Adhesive changed compared to Provox StabiliBase?

No, the finger lift tab or the pull-tab is located at the same location and has the same size as on your trusted Provox StabiliBase adhesive baseplate.

### Key features:

- ✓ Firm (yet gentle) seal and enlarged reinforced adapter for a firm good fit
- ✓ Redesigned clover-shaped adhesive for easy application and better fit
- ✓ SecureFit for security, comfort and compatibility across all Provox Life™ HMEs





QR Code, watch our online video on how to apply Provox Life™ Stability in 4 simple steps

# Top tips:

- ✓ When applying your adhesive try to fit the ring as close to the stoma as possible
- ✓ Do not forcefully pull the adhesive off as this may damage your skin
- ✓ Wait at least 5 minutes prior to speaking with a voice prosthesis to improve adhesion
- ✓ Special wipes like Provox Adhesive Remover help to peel away the adhesive easier
- ✓ Try rotating the adhesive to varying positions on the neck to find a comfortable fit for you

# Provox Life™ Sensitive

Sensitive is a gentle everyday adhesive for sensitive skin. A comfortable, gentle and low-profile adhesive that forms a good seal and provides a secure fit.



## Is Provox Life™ Sensitive Adhesive better for my skin than Provox OptiDerm?

Both Provox Life™ Sensitive and Provox OptiDerm are a skin friendly baseplate made from hypoallergenic adhesive material. Provox Life™ Sensitive is made from a new generation of hydrocolloid for a strong seal, and at the same time it is skin friendly. Provox Life™ Sensitive is also 35% thinner compared to OptiDerm to enable a better fit to different stoma shapes.

### Key features:

- ✓ With hypoallergenic material developed for sensitive skin
- ✓ Redesigned clover-shaped adhesive for easy application and better fit
- ✓ SecureFit™ for security, comfort and compatibility across all Provox Life™ HME



QR Code, watch our online video on how to apply Provox Life™ Sensitive in 4 simple steps

# Provox Life™ Adhesives for greater comfort and fit

Provox Life™ Adhesive baseplates are designed for a better fit and more comfort. They include choices suitable for a variety of skin types and stoma contours, providing you with a personalised solution.


Reasons for choosing  
Provox Life™ Adhesives:

- ✓ Softer and flexible
- ✓ Better fit
- ✓ Better adhesion



Preferred Provox Life™ Adhesives  
to their previous adhesive\*

after 12 weeks\*



## Suitable for a variety of stoma contours

Provox Life™ Adhesives are a clover-shape design which helps to fit individual stoma shapes. Each of the leaves can fold separately, which helps to provide a better fit and to contribute to an airtight adhesion.



**65%**

reported it was very easy or easy to find the right adhesive for their stoma contour\*

## Reduction in the frequency of skin irritation

Provox Life™ Sensitive and Night are designed with a hypoallergenic material which is gentle on sensitive skin.



**61%**

reported never or rarely experiencing skin irritation with Provox Life\*

## Suitable for a range of skin types

Provox Life™ Adhesive range includes three different daytime baseplates made from high-performance materials, providing you with a selection to suit your skin.



**63%**

reported it was very easy or easy to find the right adhesive for their skin type\*

## Reduction in blow-through

Provox Life™ Standard and Stability have increased glue per square meter compared with their predecessors.

This reduces the amount of mucus and air moving through the adhesive creating a more airtight seal around the stoma for speaking and everyday fit.



**36%**

Reduction in changing their adhesive due to blow-through with Provox Life\*

A simple guide to which new, improved Provox Life Adhesive might be best for you.



Provox®  
FlexiDerm™



Provox Life™  
Standard  
Adhesive



Provox®  
OptiDerm™



Provox Life™  
Sensitive  
Adhesive



Provox®  
XtraBase®



Provox®  
StabiliBase™



Provox Life™  
Stability  
Adhesive



Provox®  
Luna®  
Adhesive



Provox Life™  
Night  
Adhesive



Provox®  
StabiliBase™  
OptiDerm™



Provox Life™  
Stability  
Adhesive



Provox Life™  
Sensitive  
Adhesive





**PROVOX<sup>®</sup>**  
**Life** Breathe better,  
whatever you do

“Provox Life™ is really well thought out and great for me I have a deep irregular stoma these products fit better and reduce phlegm.”

**Provox Life™ User**

If you would like to start using Provox Life™ you can contact our Customer Care team.



**Call or text us**

Freephone: 0800 783 1659 Text: 07800 140 059



**Email us**

[info@atos-care.co.uk](mailto:info@atos-care.co.uk)



**Visit our website**

[www.atos-care.co.uk](http://www.atos-care.co.uk)



**Write to us**

Atos Care UK, Cartwright House,  
Riverside Business Park, Tottle Road,  
Nottingham, NG2 1RT

**Atos**

Breathing-Speaking-Living  
[atosmedical.com](http://atosmedical.com)

\*Provox Life Experience Program, 2021

\*\*Device specifications

